M-PRINT TryIn

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- ① en Instructions for use
- ① fr Notice d'utilisation
- ① it Istruzioni per l'uso
- ① es Instrucciones de uso
- ① pt Instruções de utilização
- ① pl Instrukcja użycia
- ① nl Gebruiksaanwijzing
- ① sv Bruksanvisning
- ① da Brugsanvisningen
- ① no Bruksanvisningen
- ① tr Kullanım kılavuzu

M-PRINT TryIn

- ① Käyttöohje
- ② Naudojimo instrukcija
- ③ Lietošanas informācija
- ④ Kasutusjuhend
- ⑤ Használati utasítás
- ⑥ Návod k použití
- ⑦ Navodila za uporabo
- ⑧ Návod na používanie
- ⑨ Инструкции за употреба
- ⑩ Οδηγίες χρήσης
- ⑪ Інструкція для застосування
- ⑫ Instructions for use

M-PRINT TryIn

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(en) Instructions for use, please read carefully!

Intended purpose

Liquid acrylate mixture for the additive fabrication of try-ins and interim prostheses in dentistry using tub-based photopolymerisation preferably at 385 nm.

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Total or partial tooth loss

Contraindication

- Hypersensitivity to one of the components
- Applications with a wear period in excess of 30 days

Composition

- Aliphatic diurethane dimethacrylate
- Ethylene dimethacrylate
- Methacrylic acid, monoester with propane-1,2-diol
- 2-[[[(butylamino)carbonyl]oxy]ethyl acrylate
- Ethyl phenyl(2,4,6-trimethylbenzoyl)phosphinate
- Phenyl bis(2,4,6-trimethylbenzoyl)phosphine oxide
- Propoxylated glycerol triacrylate
- Mequinol
- Colorants

Symbols

Color

Wavelength

Physical data

Property	Value	Test method
Acrylate mixture		
Viscosity	1200 ± 150 mPa·s	ISO 3219 (23 °C, 5 - 100 s ⁻¹)
End product		
Flexural strength.	81 ± 5 MPa	ISO 10477
Flexural strength	50 ± 2 MPa	ISO 20795-1
Modulus of elasticity	1450 ± 100 MPa	ISO 20795-1

CAD design

- Dimension the material thickness so that the wall thickness after completion is at least 1 mm for reasons of strength and a maximum of 7 mm for reasons of complete post-curing.
- When planning to bond a Bicolor construction with resin, the teeth / row of teeth should be constructed without basal cavities, as the resin present in these may not cure completely during post-curing.
- If subsequent bonding of a Bicolor construction with Weropress® is planned, select a slightly larger bonding gap than for resin bonding so that excess resin can drain well / that no increase in bite occurs and that sufficient material is available for the bond. Basal cavities for additional macro-retentions are possible.

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- A design aligned parallel to the construction platform has a lower printing time than a vertically aligned one, due to the lower number of layers, but this also requires significantly more support structures, which must be manually removed and cleaned in post-processing. It is recommended to angle the designs between 30° and 60°.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimise reworking.

Additive fabrication on a M-PRINT printer

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- The further procedure, as well as the operation of the printer, are described in the associated manual.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimise the risk of contamination and unwanted polymerisation of material still present on the surface.

Process step		Time required [min]
1	Allow the design in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (propan-2-ol; isopropanol $\geq 99\%$), such as M-PRINT Wash, for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Only perform this post-curing step if a monoblock was printed without subsequent bonding or if bonding is to be performed using Weropress® – otherwise skip this step. Post-curing of the design in the Otoflash G171 light curing device with 2 x 3000 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimises oxygen inhibition and thus leads to improved surface and material properties.	8

Bonding of Bicolor constructions

Preparation

- Cleaning of the supports with commercially available rotating instruments
- Prior to bonding, check the fit of the row of teeth in the denture base.
- The bonding surfaces of the row of teeth and denture base can be roughened using a sandblasting device and aluminium oxide (50-110 μm , 1-2 bar) to improve mechanical adhesion.
- Abrasive blasting media residues must be removed carefully using an ultrasonic bath, steam cleaner, and oil-free compressed air.

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B) Fixation with Weropress®

- Mix the Weropress® denture resin according to the IFU.
- Fill a small amount of resin into a dosing aid such as an application syringe to make it easier to apply the resin. The amount of resin depends on the size of the cavity.
- Then insert the row of teeth carefully into the corresponding cavity of the denture base until the adhesive surface is no longer visible.
- Coarse excess should be removed directly with a suitable instrument.
- This is followed by polymerisation in the pressure pot according to the IFU. Tip: using clamps, for example, to secure the dental arch, can prevent unwanted slipping of the dental arch or an increase in occlusion.

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Note: To prevent inaccurate fitting after polymerisation, avoid excessive heat development during finishing and polishing.

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Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C).
- Avoid direct sunlight.
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.

Disposal

- The fully polymerized components are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations.
Do not allow entry into the sewage system or waterways.

Notes

- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerised material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Warnings

- Processing frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the unpolymerised mixture.
- It is recommended to wear nitrile gloves when handling the liquid.
- The biocompatibility can only be assured with correct post-curing.

Hazard statements

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves and eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P273 Avoid release to the environment.

P501 Dispose of contents/container to an appropriate recycling or disposal facility.

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UA.TR.121

M-PRINT TryIn

UV-curable resin based on acrylate
for additive manufacturing

usa Instructions for use, please read carefully!

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed dental practitioner or technician.

Intended purpose

Liquid acrylate mixture for the additive fabrication of try-ins in dentistry using tub-based photopolymerization preferably at 385 nm.

Classification according to MDR (EU) 2017/745

Class I medical device for short-term contact with intact and / or injured mucosa.

Wearing time

Restricted to try-in sessions only.

The try-in device may only be worn during the patient's session at the dentist. It must not be used as a temporary, for any other purpose, or outside the designated treatment time.

Patient target groups

Adolescents and adults

Intended users

Dental professionals

Indication

Total or partial tooth loss

Contraindication

- Hypersensitivity to one of the components
- Wear duration exceeding try-in session

Composition

- Aliphatic diurethane dimethacrylate
- Ethylene dimethacrylate
- Methacrylic acid, monoester with propane-1,2-diol
- 2-[[[butylamino)carbonyl]oxy]ethyl acrylate
- Ethyl phenyl(2,4,6-trimethylbenzoyl)phosphinate
- Phenyl bis(2,4,6-trimethylbenzoyl)phosphine oxide
- Propoxylated glycerol triacrylate
- Mequinol
- Colorants

Symbols

Color

Wavelength

Physical data

Property	Value	Test method
Acrylate mixture		
Viscosity	1200 ± 150 mPa·s	ISO 3219 (23 °C, 5 - 100 s ⁻¹)
End product		
Flexural strength.	81 ± 5 MPa	ISO 10477
Flexural strength	50 ± 2 MPa	ISO 20795-1
Modulus of elasticity	1450 ± 100 MPa	ISO 20795-1

CAD design

- Dimension the material thickness so that the wall thickness after completion is at least 1 mm for reasons of strength and a maximum of 7 mm for reasons of complete post-curing.
- When planning to bond a Bicolor construction with resin, the teeth / row of teeth should be constructed without basal cavities, as the resin present in these may not cure completely during post-curing.
- If subsequent bonding of a Bicolor construction with Weropress is planned, select a slightly larger bonding gap than for resin bonding so that excess resin can drain well / that no increase in bite occurs and that sufficient material is available for the bond. Basal cavities for additional macro-retentions are possible.

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- A design aligned parallel to the construction platform has a lower printing time than a vertically aligned one, due to the lower number of layers, but this also requires significantly more support structures, which must be manually removed and cleaned in post-processing. It is recommended to angle the designs between 30° and 60°.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimize reworking.

Additive fabrication

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimize the risk of contamination and unwanted polymerization of material still present on the surface.

Process step		Time required [min]
1	Allow the design still in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (isopropanol $\geq 99\%$), for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C / 104 °F Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Only perform this post-curing step if a monoblock was printed without subsequent bonding or if bonding is to be performed using Weropress® – otherwise skip this step. Post-curing of the design in the Otoflash G171 light curing device with 2 x 3000 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimizes oxygen inhibition and thus leads to improved surface and material properties.	8

Bonding of Bicolor constructions

Preparation

- Cleaning of the supports with commercially available rotating instruments
- Prior to bonding, check the fit of the row of teeth in the denture base.
- The bonding surfaces of the row of teeth and denture base can be roughened using a sandblasting device and aluminum oxide (50 -110 μm , 1-2 bar) to improve mechanical adhesion.
- Abrasive blasting media residues must be removed carefully using an ultrasonic bath, steam cleaner, and oil-free compressed air

A) Fixation with liquid

- Fill a small amount of resin into a dosing aid such as an application syringe to make it easier to apply the resin.
- Apply the liquid to the cavity of the denture base using the dosing aid. The amount of resin depends on the size of the cavity.
- Then insert the row of teeth into the corresponding cavity of the denture base until the adhesive surface is no longer visible.
- Coarse excess should be removed directly with a suitable instrument and, if necessary, with laboratory wipes.
- This is followed by post-curing and final curing in an Otofash G171 light-curing device with 2 x 3000 flashes – rotate the construction occasionally. During the first run, align the construction such that the applied liquid cannot run out of the cavity. Tip: using clamps, for example, to secure the dental arch, can prevent unwanted slipping of the dental arch or an increase in occlusion.

B) Fixation with Weropress

- Mix the Weropress denture resin according to the IFU.
- Fill a small amount of resin into a dosing aid such as an application syringe to make it easier to apply the resin. The amount of resin depends on the size of the cavity.
- Then insert the row of teeth carefully into the corresponding cavity of the denture base until the adhesive surface is no longer visible.
- Coarse excess should be removed directly with a suitable instrument.
- This is followed by polymerization in the pressure pot according to the IFU. Tip: using clamps, for example, to secure the dental arch, can prevent unwanted slipping of the dental arch or an increase in occlusion.

Manual finishing & polishing

Manual post-processing, e.g. cleaning of the supports and polishing with commercially available rotating instruments for plastic processing.

Note: To prevent inaccurate fitting after polymerization, avoid excessive heat development during finishing and polishing.

Disinfection

- The post-cured designs can be disinfected with both alcoholic and non-alcoholic disinfectants.
- Do not use heat-based methods, as this can lead to deformations.

Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C / 59 °F to 77 °F).
- Avoid direct sunlight.
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.

Disposal

- The fully polymerized components are disposed of as residual waste.
- Dispose of non-polymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerized material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Warnings

- Processing frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the unpolymerized mixture.
- It is recommended to wear nitrile gloves when handling the liquid.
- The biocompatibility can only be assured with correct post-curing.

Hazard statements

Causes skin irritation and serious eye irritation

May cause an allergic skin reaction

Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Contaminated work clothing must not be allowed out of the workplace.

Wear protective gloves and eye protection/face protection.

If on skin: Wash with plenty of water.

If skin irritation or rash occurs: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists: Get medical advice/attention.

Dispose of contents/container to an appropriate recycling or disposal facility.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from **www.merz-dental.de**.

The product properties are based on compliance with and observation of these instructions for use.

The users themselves are responsible for the use of the product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

Date of information 2025-05

Merz Dental GmbH

Kieferweg 1, 24321 Lütjenburg, Germany

Tel + 49 (0) 4381 / 403-0

Fax + 49 (0) 4381 / 403-403

www.merz-dental.de

EN ISO 13485

907402/2025-05/01