M-PRINT Model 2.0

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USA **en** Instructions for use

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M-PRINT Model 2.0

UV-curable resin based on acrylate
for additive manufacturing

en Instructions for use, please read carefully!

Intended purpose

Liquid acrylate mixture (Liquid) for the additive manufacturing of models in dentistry (final product) using tub-based photopolymerisation at 385 - 405 nm.

Composition of the liquid	desert & shark-gray	deep-black
Diurethane dimethacrylate, mix of isomers (UDMA)	X	X
(3-Phenoxyphenyl)methyl 2-propenoate	X	X
Silicone Surfactants	X	-
Ethylene dimethacrylate	X	-
2-[[[(butylamino)carbonyl]oxy]ethyl acrylate	X	X
Mequinol	X	X
Colorants	X	X

Symbols

Color

Wavelength

Physical data

Property	Value		Test method
Acrylate mixture	desert & shark-gray	deep-black	
Viscosity	1250 ± 150 mPa·s	1250 ± 150 mPa·s	ISO 3219
End product			
Flexural strength	90 ± 5 MPa	128 ± 5 MPa	ISO 178
Modulus of elasticity	2150 ± 150 MPa	3200 ± 200 MPa	ISO 178

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 25 µm, 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimise reworking.

Additive fabrication on a M-PRINT printer

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- The further procedure, as well as the operation of the printer, are described in the associated manual.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimise the risk of contamination and unwanted polymerisation of material still present on the surface.

Process step		Time required [min]
1	Allow the design still in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (propan-2-ol; isopropanol $\geq 99\%$), such as M-PRINT Wash, for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Post-curing of the design in the Otoflash G171 light curing device with 2 x 3500 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimises oxygen inhibition and thus leads to improved surface and material properties. Observe the manufacturer's specifications for the device.	8
6	Manual post-processing, e.g. cleaning of the supports and polishing with commercially available rotating instruments for plastic processing. Note: To prevent inaccurate fitting after polymerisation, avoid excessive heat development during finishing and polishing.	Design-dependent

Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C).
- Avoid direct sunlight.
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.
- Long-term storage of the liquid in the material tub can reduce its service life.
- For short-term storage, such as when changing materials, it is advisable to store the material-filled tub in the original packaging, as this protects against dirt and light at the same time.

Disposal

- The fully polymerized final product or its components are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Do not use the liquid after the expiration date.
- Record the LOT number with every process that requires identification of the liquid.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerised material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.

Warnings

- Processing final products may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the liquid.
- It is required to wear nitrile gloves when handling the liquid.

Hazard statements (Liquid)

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H412 Harmful to aquatic life with long lasting effects. (deep-black)

Precautionary statements (Liquid)

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing and eye protection/face protection.

P362+P364 Take off contaminated clothing and wash it before reuse.

P501 Dispose of contents/container to an appropriate recycling or disposal facility.

For detailed information about the liquid please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with and observation of these instructions for use.

The users themselves are responsible for the proper use of the liquid and the use of the final product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

Date of information 2026-01

M-PRINT Model 2.0

UV-curable resin based on acrylate
for additive manufacturing

(USA) Instructions for use, please read carefully!

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed dental practitioner and dental technician.

Intended purpose

Liquid acrylate mixture (Liquid) for the additive manufacturing of models in dentistry (final product) using tub-based photopolymerisation at 385 - 405 nm.

Composition of the liquid	desert & shark-gray	deep-black
Diurethane dimethacrylate, mix of isomers (UDMA)	X	X
(3-Phenoxyphenyl)methyl 2-propenoate	X	X
Silicone Surfactants	X	-
Ethylene dimethacrylate	X	-
2-[[[(butylamino)carbonyl]oxy]ethyl acrylate	X	X
Mequinol	X	X
Colorants	X	X

Symbols

Color

Wavelength

Physical data

Property	Value		Test method
Acrylate mixture	desert & shark-gray	deep-black	
Viscosity	1250 ± 150 mPa·s	1250 ± 150 mPa·s	ISO 3219
End product			
Flexural strength	90 ± 5 MPa	128 ± 5 MPa	ISO 178
Modulus of elasticity	2150 ± 150 Mpa	3200 ± 200 Mpa	ISO 178

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 25 µm, 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimise reworking.

Additive fabrication on a M-PRINT printer

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- The further procedure, as well as the operation of the printer, are described in the associated manual.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimise the risk of contamination and unwanted polymerisation of material still present on the surface.

Process step		Time required [min]
1	Allow the design still in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (propan-2-ol; isopropanol $\geq 99\%$) for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C / 104 °F Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Post-curing of the design in the Otofash G171 light curing device with 2 x 3500 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimises oxygen inhibition and thus leads to improved surface and material properties. Observe the manufacturer's specifications for the device.	8
6	Manual post-processing, e.g. cleaning of the supports and polishing with commercially available rotating instruments for plastic processing. Note: To prevent inaccurate fitting after polymerisation, avoid excessive heat development during finishing and polishing.	Design-dependent

Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C / 59 °F to 77 °F).
- Avoid direct sunlight.
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.
- Long-term storage of the liquid in the material tub can reduce its service life.
- For short-term storage, such as when changing materials, it is advisable to store the material-filled tub in the original packaging, as this protects against dirt and light at the same time.

Disposal

- The fully polymerized final product or its components are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Do not use the liquid after the expiration date.
- Record the LOT number with every process that requires identification of the liquid.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerised material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.

Warnings

- Processing final products may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the liquid.
- It is required to wear nitrile gloves when handling the liquid.

Hazard statements (Liquid)

Causes skin irritation and serious eye irritation.

May cause an allergic skin reaction.

Precautionary statements (Liquid)

Avoid breathing dust/fume/gas/mist/vapors/spray.

Contaminated work clothing must not be allowed out of the workplace.

Wear protective gloves/protective clothing and eye protection/face protection.

If skin irritation or rash occurs: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

Dispose of contents/container to an appropriate recycling or disposal facility.

For detailed information about the liquid please see the applicable Material Safety Data Sheet, available upon request or as download from **www.merz-dental.de**.

The product properties are based on compliance with and observation of these instructions for use.

The users themselves are responsible for the proper use of the liquid and the use of the final product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

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Merz Dental GmbH