M-PRINT

Model ISO

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- | | |
|-------------------------------|----------------------------------|
| (en) Instructions for use | (lt) Naudojimo instrukcija |
| (fr) Notice d'utilisation | (lv) Lietošanas informācija |
| (it) Istruzioni per l'uso | (et) Kasutusjuhend |
| (es) Instrucciones de uso | (hu) Használati utasítás |
| (pt) Instruções de utilização | (cs) Návod k použití |
| (pl) Instrukcja użytkowania | (sl) Navodila za uporabo |
| (nl) Gebruiksaanwijzing | (sk) Návod na používanie |
| (sv) Bruksanvisning | (bg) Инструкции за употреба |
| (da) Brugsanvisningen | (el) Οδηγίες χρήσης |
| (no) Bruksanvisningen | (tr) Kullanım kılavuzu |
| (fi) Käyttöohje | (uk) Інструкція для застосування |

M-PRINT Model ISO

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M-PRINT Model ISO

Ready-to-use isolating agent for 3D printed resin-models

(en) Instructions for use, please read carefully!

Intended purpose

Ready-to-use isolating agent for 3D-printed resins - e.g. against denture resins (cold-curing resin) as well as composite veneering materials.

Composition

Propan-2-ol

Preparation

- 3D Resins must be fully post cured, cleaned, and dry (Incomplete polymerized 3D resins may impair the insulating effect).
- Before application of the isolating agent, block out undercuts with wax. The isolating agent also sticks to the wax.

Application

- Apply M-PRINT Model ISO once, using either a natural or synthetic brush.
- Coat M-Print Model ISO directly onto the model and apply it a bit thicker (especially in filament models), afterwards the isolation can be coated into a uniform layer. This procedure does not leave any uninsulated areas, the isolation agent does not dry as quickly when applied and prevents the dissolution of dried layers. Excess material can be easily and accurately blown off with oil-free compressed air or wiped off with a paper towel.
- Disperse any puddles (e.g., in deep fissures).
- Leave it to dry completely.

Cleaning:

- Rinse the brush with clear water
- The isolation layer can be removed with water or by steamcleaning.

Storage

Shelf life at least 24 months in closed containers.

Disposal

Dispose of contents/containers in accordance to local regulations.

Instructions and tips

- Keep the bottle closed tightly after use.
- To avoid contamination and depositing dry isolation, we recommend applying directly onto the model.

- Applying a second layer is not necessary for models produced using the DLP process. In the case of a second application, the dried isolation layer may become brittle and lead to a deterioration of the surface (Chequered appearance).
- For models fabricated from filament or in case of pronounced layers, repeated application may be helpful on account of the surface quality. Hardening in a light curing device is recommended for improved drying.
- The blue colour scheme allows a good visual inspection of the isolated areas and the layer thickness. The dried isolating area is easily recognizable.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

Hazard statements

H225 Highly flammable liquid and vapour.

H319 Causes serious eye irritation.

H336 May cause drowsiness or dizziness.

Precautionary statements

P233 Keep container tightly closed.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

The product properties are based on compliance with and observation of these instructions for use.

Printed copies of these instructions for use can be requested free of charge using the contact details provided.

The users themselves are responsible for the use of the product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

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Merz Dental GmbH