M-PRINT Splint

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- ① **USA** ② **en** Instructions for use
- ③ **fr** Notice d'utilisation
- ④ **it** Istruzioni per l'uso
- ⑤ **es** Instrucciones de uso
- ⑥ **pt** Instruções de utilização
- ⑦ **nl** Gebruiksaanwijzing
- ⑧ **el** Οδηγίες χρήσης
- ⑨ **tr** Kullanım kılavuzu

M-PRINT Splint

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M-PRINT Splint

UV-curable resin based on acrylate
for additive manufacturing

en Instructions for use, please read carefully!

Intended purpose

Liquid acrylate mixture (Liquid) for the additive fabrication of functional and impression trays, as well as bite templates in dentistry (final product) using tub-based photopolymerisation, preferably at 385 nm.

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Liquid acrylate mixture (Liquid): Medical device as an intermediate product of class IIa.

Final product to be manufactured: Class IIa medical device for long-term contact with intact mucous membrane.

Patient target groups

Adolescents and adults

Intended users

Dental professionals

Indication of the final product

Total or partial tooth loss

Contraindication of the final product

Hypersensitivity to one of the components

Composition of the liquid

- Diurethane dimethacrylate (UDMA)
- 2-[[[(butylamino)carbonyl]oxy]ethyl]acrylate
- 1,4-Butanediol dimethacrylate
- Urethane acrylate in polyether polyol tetraacrylate
- Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide (TPO)
- Propoxylated glycerol triacrylate
- Mequinol
- Colorants

Symbols

Color

Wavelength

Physical data

Property	Value	Test method
Acrylate mixture		
Viscosity	825 ± 100 mPa·s	ISO 3219 (23 °C, 5 - 100 s ⁻¹)
End product		
Flexural strength	70 ± 6 MPa	ISO 20795-1
Modulus of elasticity	1800 ± 200 MPa	ISO 20795-1
Glass transition	93 °C	ISO 11357
Vickers hardness	20 ± 1 HV0,2	ISO 6507-1*

*based on

CAD design

Dimension the material thickness so that the wall thickness after completion is at least 1,5 mm for reasons of strength and a maximum of 7 mm for reasons of complete postcuring.

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 25 µm, 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- A design aligned parallel to the construction platform has a lower printing time than a vertically aligned one, due to the lower number of layers, but this also requires significantly more support structures, which must be manually removed and cleaned in post-processing. It is recommended to angle the designs between 30° and 90°.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimise reworking.

Additive fabrication on a M-PRINT printer

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- The further procedure, as well as the operation of the printer, are described in the associated manual.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimise the risk of contamination and unwanted polymerisation of material still present on the surface.

Process step		Time required [min]
1	Allow the design still in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (propan-2-ol; isopropanol $\geq 99\%$), such as M-PRINT Wash, for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Post-curing of the design in the Otofash G171 light curing device with 2 x 3000 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimises oxygen inhibition and thus leads to improved surface and material properties. Observe the manufacturer's specifications for the device.	8
6	Manual post-processing, e.g. cleaning of the supports and polishing with commercially available rotating instruments for plastic processing. Note: To prevent inaccurate fitting after polymerisation, avoid excessive heat development during finishing and polishing.	Design-dependent

Disinfection of the final product

- The post-cured final product can be disinfected with both alcoholic and non-alcoholic disinfectants.
- Do not use heat-based methods, as this can lead to deformations.

Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C).
- Avoid direct sunlight.
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.
- Long-term storage of the liquid in the material tub can reduce its service life.
- For short-term storage, such as when changing materials, it is advisable to store the material-filled tub in the original packaging, as this protects against dirt and light at the same time.

Disposal

- The fully polymerized final product or its components are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Do not use the liquid after the expiration date.
- Record the LOT number with every process that requires identification of the liquid.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerised material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.
- All serious incidents occurring with users and/or patients in connection with the liquid or the final product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and/or patient is resident.

Warnings

- Processing final products may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the liquid.
- It is required to wear nitrile gloves when handling the liquid.
- The biocompatibility of the final product can only be assured with correct post-curing.

Hazard statements (Liquid)

H317 May cause an allergic skin reaction.

H360F May damage fertility.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements (Liquid)

P202 Do not handle until all safety precautions have been read and understood.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing and eye protection/face protection.

P362+P364 Take off contaminated clothing and wash it before reuse.

P273 Avoid release to the environment.

P501 Dispose of contents/container to hazardous or special waste collection point.

For detailed information about the liquid please see the applicable Material Safety Data Sheet, available upon request or as download from **www.merz-dental.de**.

The product properties are based on compliance with and observation of these instructions for use.

The users themselves are responsible for the proper use of the liquid and the use of the final product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

Date of information 2025-09

M-PRINT Splint

UV-curable resin based on acrylate
for additive manufacturing

USA Instructions for use, please read carefully!

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed dental practitioner and dental technician.

Intended purpose

Liquid acrylate mixture (Liquid) for the additive fabrication of functional and impression trays, as well as bite templates in dentistry (final product) using tub-based photopolymerisation, preferably at 385 nm.

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Liquid acrylate mixture (Liquid): Medical device as an intermediate product of class IIa.

Final product to be manufactured: Class IIa medical device for long-term contact with intact mucous membrane.

Patient target groups

Adolescents and adults

Intended users

Dental professionals

Indication of the final product

Total or partial tooth loss

Contraindication of the final product

Hypersensitivity to one of the components

Composition of the liquid

- Diurethane dimethacrylate (UDMA)
- 2-[[[(butylamino)carbonyl]oxy]ethyl]acrylate
- 1,4-Butanediol dimethacrylate
- Urethane acrylate in polyether polyol tetraacrylate
- Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide (TPO)
- Propoxylated glycerol triacrylate
- Mequinol
- Colorants

Symbols

Color

Wavelength

Physical data

Property	Value	Test method
Acrylate mixture		
Viscosity	825 ± 100 mPa·s	ISO 3219 (23 °C, 5 - 100 s ⁻¹)
End product		
Flexural strength	70 ± 6 MPa	ISO 20795-1
Modulus of elasticity	1800 ± 200 MPa	ISO 20795-1
Glass transition	93 °C / 199,4 °F	ISO 11357
Vickers hardness	20 ± 1 HV0,2	ISO 6507-1*

*based on

CAD design

Dimension the material thickness so that the wall thickness after completion is at least 1,5 mm for reasons of strength and a maximum of 7 mm for reasons of complete postcuring.

CAM processing

- Use the correct printer / material combination. Released parameter sets are available for download at www.merz-dental.de
- The material is designed for layer thicknesses / slices of 25 µm, 50 µm and 100 µm.
- To avoid dimensional errors in the z-axis, it is recommended to use support structures to place the design on a base plate with a perforated pattern, e.g. a hexagonal grid.
- A design aligned parallel to the construction platform has a lower printing time than a vertically aligned one, due to the lower number of layers, but this also requires significantly more support structures, which must be manually removed and cleaned in post-processing. It is recommended to angle the designs between 30° and 90°.
- The use of predetermined breaking points / tapering of the supports towards the component is recommended, as this reduces the risk of damaging the surface when removing the supports.
- When positioning the supports, make sure that overhangs are supported and, wherever possible, fitting surfaces are not pinned to minimise reworking.

Additive fabrication on a M-PRINT printer

- Use the parameter set corresponding to the material.
- For printing, use the same layer thickness as in CAM planning, otherwise inaccuracies in fit and / or misprints may occur.
- Ensure clean work. Contamination, especially at the optical window to the projector and at the material tub (vat), can cause errors on the printed part or misprints.
- When filling the material tub, make sure that the material is exposed to ambient light as briefly as possible and that the filling level markings are observed.
- The further procedure, as well as the operation of the printer, are described in the associated manual.
- After printing, the design should remain suspended in the printer for 10 minutes to allow excess print material to drip off the surface. This saves material and prolongs the service life of the cleaning fluid.

Post-processing

Post-processing should be performed immediately after completion of printing to minimise the risk of contamination and unwanted polymerisation of material still present on the surface.

Process step		Time required [min]
1	Allow the design still in the printer after the printing process has been completed to drip off.	10
2	Remove construction platform, remove design and carefully cut off supports.	5
3	Clean the design in a closed container with cleaning fluid (propan-2-ol; isopropanol $\geq 99\%$) for 5 minutes in an ultrasonic bath. Note: Do not heat the design in the cleaning solution and do not leave it in the solution for a long period.	5
4	Use a pair of tweezers to remove the printed part from the cleaning fluid, blow it off with compressed air if necessary, and let it dry. Recommendation: Dry the design in an oven at 40 °C / 104 °F Note: Residues of the alcoholic cleaning fluid on the surface lead to a softer and sometimes scratch-sensitive surface after post-curing.	30
5	Post-curing of the design in the Otofash G171 light curing device with 2 x 3000 flashes - turn the design between the two cycles. Recommendation: Use the nitrogen (protective gas) function, as this minimises oxygen inhibition and thus leads to improved surface and material properties. Observe the manufacturer's specifications for the device.	8
6	Manual post-processing, e.g. cleaning of the supports and polishing with commercially available rotating instruments for plastic processing. Note: To prevent inaccurate fitting after polymerisation, avoid excessive heat development during finishing and polishing.	Design-dependent

Disinfection of the final product

- The post-cured final product can be disinfected with both alcoholic and non-alcoholic disinfectants.
- Do not use heat-based methods, as this can lead to deformations.

Storage of the liquid

- Store dry and at room temperature (15 °C to 25 °C / 59 °F to 77 °F).
- Avoid direct sunlight
- Keep containers closed and seal carefully after use.
- Do not return residual amounts from the material tub to the original container.
- Long-term storage of the liquid in the material tub can reduce its service life.
- For short-term storage, such as when changing materials, it is advisable to store the material-filled tub in the original packaging, as this protects against dirt and light at the same time.

Disposal

- The fully polymerized final product or its components are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Do not use the liquid after the expiration date.
- Record the LOT number with every process that requires identification of the liquid.
- When selecting the parameter set, make sure it matches the version number of the liquid on the batch label.
- Due to the reactivity of the liquid, it is recommended to clean wetted parts, such as forceps and the construction platform, immediately after use, as the polymerised material is difficult to remove.
- If the cleaning fluid becomes noticeably cloudy, the cleaning capacity is spent and it should be replaced with fresh fluid. If there are still damp, sticky patches on the designs after the drying process, this is also an indication of a saturated cleaning fluid and the design should be cleaned again with fresh cleaning fluid.
- When changing material / cleaning the material tub, use soft cloths and the silicone wiper M-PRINT Wipe to avoid damaging the surface of the material tub.
- If other printers are used, make sure that the appropriate printer data sets are available in the CAM and that material parameters for the printers are available.
- All serious incidents occurring with users and/or patients in connection with the liquid or the final product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and/or patient is resident.

Warnings

- Processing final products may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- Avoid direct skin contact with the liquid.
- It is required to wear nitrile gloves when handling the liquid.
- The biocompatibility of the final product can only be assured with correct post-curing.

Hazard statements (Liquid)

May cause an allergic skin reaction

May damage fertility; suspected of damaging the unborn child

Precautionary statements (Liquid)

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

Avoid breathing dust/fume/gas/mist/vapors/spray.

Contaminated work clothing must not be allowed out of the workplace.

Wear protective gloves/protective clothing and eye protection/face protection.

If on skin: Wash with plenty of water.

Wash contaminated clothing before reuse.

If exposed or concerned: Get medical advice/attention.

Store locked up.

Dispose of contents/container to hazardous or special waste collection point.

For detailed information about the liquid please see the applicable Material Safety Data Sheet, available upon request or as download from **www.merz-dental.de**.

The product properties are based on compliance with and observation of these instructions for use.

The users themselves are responsible for the proper use of the liquid and the use of the final product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

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Merz Dental GmbH