PremEco® Line

Alginatisolierung

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- | | |
|--------------------------------------|------------------------------------|
| (en) Instructions for use | (fi) Käyttöohje |
| (fr) Notice d'utilisation | (it) Naudojimo instrukcija |
| (it) Istruzioni per l'uso | (lv) Lietošanas informācija |
| (es) Instrucciones de uso | (et) Kasutusjuhend |
| (pt) Instruções de utilização | (hu) Használati utasítás |
| (pl) Instrukcja użycia | (cs) Návod k použití |
| (nl) Gebruiksaanwijzing | (sl) Navodila za uporabo |
| (sv) Bruksanvisning | (sk) Návod na používanie |
| (da) Brugsanvisningen | (el) Οδηγίες χρήσης |
| (no) Bruksanvisningen | (tr) Kullanım kılavuzu |

PremEco® Line

Alginatisolierung

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Stand der Information 2013-05

PremEco® Line Alginate isolating solution

Special isolating solution to use when processing heat-curing and cold-curing polymers. Thanks to its optimum consistency it can be easily applied, absorbs easily into plaster and isolates well.

(en) Instructions for use, please read carefully!

Areas of Application

Isolation for all completion work, relines, repairs and additions to partial and full dentures, saddles of model cast dentures, removable saddles of combined dental restorations, immediate and temporary dentures.

Composition

PremEco® Line Alginate isolating solution is a formaldehyde-free mixture of water, alginate and preservatives.

Dosage and Method of Administration

- Clean off, separating agent and grease from models and moulds and wet thoroughly with water.
- Allow the plaster surface to dry briefly.
- After boiling out allow the models to cool until water vapour is no longer released.
- Fill the required amount of the insulator into a container.
- Apply the alginate isolating solution once or twice thinly using a clean, soft brush.
- Never dip the brush directly into the bottle containing the alginate insulation solution.
- In order to prevent contamination of the material, do not return any residual quantities back into the bottle.
- Clean used brushes under running water.
- Do not apply to clasps framework parts and teeth.
- Avoid excess solution in grooves, furrows and folds.
- Rinse the isolated model carefully under running water.
- Satisfactory isolation can be identified by a uniform gloss on the surface.

PremEco® Line Alginate isolating solution is a product from the PremEco® Line fabrication system.

For further information please visit our website at www.merz-dental.de.

Special Notes

- To prevent contamination, carefully close the PremEco® Alginate isolating solution bottle after use.
- Store in a dry place at ambient temperature; avoid direct exposure to sunlight.
- Do not use the product after the expiry date.
- Indicate the LOT No. in all processes requiring identification of the material.

The product properties are based on compliance with and observation of these instructions for use.

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Merz Dental GmbH