UBIGEL Inperio

CE 0476 - Dispositivo Medico

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UBIGEL Inperio

CE 0476 - Medical Device

Periodontal gel - Periodontal regenerating adhesive as an adjuvant in the treatment of gingivitis and/or periodontitis

UBIGEL Inperio is a periodontal regenerating adhesive gel indicated as an adjuvant in the treatment of gingivitis and/or periodontitis.

INSTRUCTIONS: *UBIGEL Inperio* is a filling gel for periodontal pockets in gums affected by periodontitis. It improves and supports the local trophism of tissues affected by gingivitis and periodontitis, promoting their reparative processes.

Use is reserved for dentists. Periodontitis is a disease in which the action of some bacteria leads to the inflammatory destruction of the supporting tissues of the teeth, with consequent recession of the gums and resorption of the bone. The new therapies promote the healing of periodontitis by arresting its progression and sometimes allowing the regeneration of periodontal tissues. *UBIGEL Inperio* is a gel containing an addition complex with spermidine, alginate, and hyaluronic acid with a trophic action. Alginate has a dual function: a powerful adhesive agent and a controlled release matrix of spermidine to optimise its trophic activity, creating the optimal condition for restoring the physiological conditions of the oral mucosa. Spermidine is a polyamine with a ubiquitous character present in the tissues, with a high capacity for regeneration of the connective, bone and dental tissue damaged by the erosive effects of periodontitis. It therefore promotes cell growth and differentiation, useful in the repair of mucous membranes. The *Ubigel Inperio* medical device is capable of creating a "soft cement" within the periodontal pocket, of the consistency of a rubber, with a certain solidity to perform the following physical actions:

1. occlusive action with the aim of preventing the infiltration of food residues that increase the infection;
2. mechanical abrasive action on the bacterial film of the periodontal pocket to reduce or eliminate its presence;
3. sealing action against bacterial spread thanks to the muco-adhesive properties of the polymeric components;
4. antibacterial action.

The nature of *UBIGEL Inperio* means it can be used in conjunction with treatments aimed at oral infections.

COMPOSITION: Spermidine; Sodium Alginate; Sodium hyaluronate; Chlorobutanol; Calcium chloride; Xylitol; Benzyl alcohol; Phenoxyethanol;

INSTRUCTIONS FOR USE: Open the package and take out the two syringes and the 2 sterile irrigation needles with side hole.

Step one: take the syringe with trophic action gel A, BLUE colour, shake vigorously, unscrew the cap and fit one of the two needles, insert the needle into the periodontal pocket until it touches the bottom and dispense the gel until the pocket is full.

Step two: take the other syringe with sealing gel B, RED colour, unscrew the cap and mount the remaining needle, dispense enough gel to close the periodontal pocket already filled with Gel A.

Check the procedure and repeat the infiltration if no adhesion. For the treatment to work more effectively, and to ensure adhesion of the gel in the intervention area, we recommend adequate scaling and curettage procedures of the periodontal pocket(s) before using *Ubigel Inperio*.

Caution: The needles in the package are EO sterile (treatment with ethylene oxide), for correct use follow the instructions given in the pictograms. If the seal is broken, the needle is NOT STERILE and MUST NOT BE USED. This device is intended for single patient use. Any re-use exposes other patients to critical risks, including possible serious infections.

FREQUENCY OF USE AND SUGGESTED DOSES: Proceed to infiltrate the gel into the periodontal pockets affected by periodontitis, from 1 to a maximum of 3 infiltrations, each one week apart. The 0.25 mL strength (single-dose, single use) corresponds to the amount needed to fill a periodontal pocket. The 1 mL dosage (multi-dose, single use) is for the treatment of multiple gingival sites of the same patient in the same treatment session.

WARNINGS: Postpone oral hygiene for 12 hours after treatment with *UBIGEL Inperio* and avoid using dental floss in the treated area for about 2 weeks. Advise the patient to avoid hard, crunchy or sticky foods that can cause trauma to the gum. Do not use the medical device after the expiry date shown on the package. Keep the medical device away from heat sources. Do not use the device if the external packaging or the security seal are damaged upon first use.

INTERACTIONS: There are currently no known interactions between *UBIGEL Inperio* and medicinal products normally used in relevant drug therapies. The product is compatible with normal toothpaste and other oral hygiene products.

CONTRAINDICATIONS AND SIDE EFFECTS: Known individual hypersensitivity to the components. Since complete data is not available, the product is not recommended for children except in cases of real need and under close medical supervision.

VALIDITY: Preferably use within 24 months

STORAGE: Keep the product at room temperature, storage in the refrigerator at 4-8 °C is the best for product conservation.

UBIGEL Inperio PACKAGING: The pack can be 0.25 mL or 1.00 mL

Box containing 2 syringes (1 Gel A BLU colour + 1 Gel B RED colour) each pre-filled with 0.25 mL of gel (single use - single-dose)

Box containing 2 syringes (1 Gel A BLU colour + 1 Gel B RED colour) each pre-filled with 1.00 mL of gel (single use - multi-dose)

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