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(en) Instructions for use, please read carefully!

Regulatory information

^{mp}Illusion Color and ^{mp}Illusion Paste meet all requirements of applicable directives and regulations on medical devices, and the manufacturing is certified acc. ISO 13485 and annex II of MDD 93/42/EWG, VO (EU) 2017/745 and further international requirements.

Medical device classification acc. annex IX, rule 8 of MDR 93/42: IIa

Medical device classification acc. annex VIII, rule 8 of MDR 2017/745: IIa

MDR Code acc. MDCG 2019-14: MDT 2003, MDN 1103

Classification acc. DIN EN ISO 6872: type 1, class 1

Intended use in general:

Dental ceramics for the fabrication of all-ceramic (crowns, bridges, inlays, onlays, veneers) or metal-ceramic restorations (crowns, bridges).

Product description

The pasty ^{mp}Illusion Color and ^{mp}Illusion Paste are used for individual characterization in the fabrication of esthetically high-quality dental restorations.

Scope

Reconstruction of natural teeth (crowns and bridges) made from following framework materials:

1. monolithic or layered zirconium dioxide (Y-TZP) crowns and bridges with thermal expansion of $10.6 \times 10^{-6} \times K^{-1}$ (25-500 °C).
2. Lithium disilicate (monolithic / veneered) with a thermal expansion of approx. $10.0 \times 10^{-6} \times K^{-1}$ (25-500 °C).
3. Frameworks made of pressed ceramics with a CTE value of approx. $9.5 \times 10^{-6} \times K^{-1}$
4. Metal ceramic with a CTE value of approx. 13.0 ± 0.5 (25-500 °C) [$10^{-6} K^{-1}$].

Patient target group

Patients with retained teeth

Intended users/special training

1. Dentist (clinical work procedure)
 2. Dental technician (fabrication of the restoration)
- No special training required

Use

Exclusively by dental professionals

Restrictions on use

1. Combinations with materials outside of the product system described.
2. Veneering of framework materials that are not indicated.
3. Sharp edges and corners on the framework.
4. Dental ceramic and all ceramic restorations made of glass ceramics are not recommended for patients with bruxism or parafunction.
5. The final restoration must not be reused.

Chemical composition ^{mp}Illusion Paste:

SiO₂, Al₂O₃, K₂O, Na₂O, Li₂O, SrO, B₂O₃, CeO₂, ZnO + pigments

Chemical composition ^{mp}Illusion Color:

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Types of restorations:

1. Crowns
2. Bridges
3. Veneers, inlays, onlays

Contraindications

For proven allergies to the ingredients contained in this product.

Adverse reactions

No adverse reactions are known at this time.

Interactions

No interactions are known at this time.

Residual risks

Users should be aware that dental interventions in the oral cavity generally involve certain risks. Some examples are given in the following:

1. Chipping/fracture/decementation of the restorative material can lead to the swallowing or inhaling of material and the need for further dental treatment.
2. Excess cement can irritate the tissue/gingiva. Continued inflammation can lead to bone loss and periodontitis.

Explanation of symbols on the product labels

Manufacturer

Batch code

Date of manufacture YYYY-MM

Reference number

Medical device Rx Only

eIFU

Caution, consult instruction for use

www.merz-dental.de/eifu

Physical properties

^{my}Illusion Paste

Physico-chemical properties according to DIN EN ISO 6872

Property	Specification	Measurement values
Coefficient of thermal expansion (25 - 475 °C) [$\cdot 10^{-6} \cdot K^{-1} \pm 0.5$]	2 x: 9.5 4 x: 9.5	2 x: 9.5 4 x: 9.5
Glass transition temperature T _g [°C ± 20]	2 x: 495 4 x: 495	2 x: 495 4 x: 495
Flexural strength [MPa]	≥ 50	130 - 150
Solubility [$\mu g/cm^2$]	< 100	19 - 35

^mIllusion Color

Physico-chemical properties according to DIN EN ISO 6872

Property	Specification	Measurement values
Coefficient of thermal expansion (25 - 450 °C) [$10^{-6} \cdot K^{-1} \pm 0.5$]	2 x: 10.0 4 x: 10.0	2 x: 9.8* 4 x: 9.5*
Glass transition temperature T _g [°C ± 20]	2 x: 460 4 x: 460	2 x: 455* 4 x: 455*
Flexural strength [MPa]	≥ 50	> 130*
Solubility [µg/cm ²]	< 100	complies with*

*) data for base material

Important processing instructions

- Only to be used by trained personnel.
- For use in clean working environments only! Contamination of the pastes, the mixing tray, the preheating furnace or any additional materials as waxes or liquids especially with CoCr-alloy residues may cause discoloration of restorations.

Recommended storage conditions

- 12 - 38° C and normal air humidity 40 - 60%.
- Storage in tightly closed original containers.
- Protect from direct sunlight.

Application instructions

- Before each application, mix the ^mIllusion Color and ^mIllusion Paste compounds well with a metal-free spatula.
- Always use a clean and dry spoon, spatula or brush.
- Moisten brush with ^mIllusion Paste liquid before application.
- Only use ^mIllusion Paste Liquid to adjust the consistency of both the ^mIllusion Color and ^mIllusion Paste compounds.
- **Note:** The pastes must not come into contact with water.
- Do not fill removed and used pastes back into the container. Use clean and dry instruments for removal.

Warnings

- When working on ceramic restorations safety glasses should be used. Remove dust and fragments by suction or wear a proper dust protection mask.
- Be careful of high firing and pressing temperatures. Danger of getting burnt! Use oven pincers and gloves!

Framework preparation

- Prepare the frameworks according to the manufacturer's instructions.
- Before the first application of ^mIllusion Color or ^mIllusion Paste, carefully sandblast the surface with 110 µm aluminum oxide and 2 bar pressure.
- Steam clean the objects before each application.

Colouring

- Before colouring, apply a thin layer of ^mIllusion Liquid to the surface.

Structuring

- Before applying ^mIllusion Paste for the first time, apply a thin layer of ^mIllusion Paste Liquid to the surface.

Layering/staining

Apply ^mIllusion product components to the crown or bridge as desired.

Applying the glaze/polishing

After finishing the restoration, apply a thin layer of ^mIllusion glaze fluor paste to the entire crown and fire. The gloss level of the restoration can be adjusted after the glaze firing using a polishing wheel and polishing paste.

FIRING TABLE

Due to the different ceramic furnaces available on the market, the firing conditions may differ. This must be taken into account and is under the responsibility of the customer!

The specified firing temperatures are **only reference values!**

Note: The indicated firing temperatures are only APPROXIMATE VALUES tested in a Zubler Vario 300 ceramic furnace. Adjustments to the firing program may be necessary for other furnace types.

Firing parameters	Start-temperature	Pre-Dry time	Closing time	Heating rate	Vakuum start	Final temperature	Holding time
	(°C)	(min)	(min)	(°C/min)	(°C)	(°C)	(min)
Monolithic zirconium dioxide (Y-TZP)	400	3	4	45	670	770	1
Metal ceramic	400	3	4	45	670	770	1
Pressed ceramic	400	3	4	45	670	770	1
Lithium disilicate (monolithic)	400	3	4	45	670	720	2
Lithium disilicate (veneered)	400	3	4	45	670	720	2
Zirconium dioxide (Y-TZP), veneered	400	3	4	45	670	720	2

Please note: please use retarded furnace opening time (2 min.) for voluminous reconstructions!

Safety notes

In the event of serious events which have occurred in connection with the product, please contact estetic ceram ag, Landstrasse 109, 9495 Triesen, Liechtenstein and your local state authority.

A brief report on the Summary of Safety and Clinical Performance (SSCP) can be found on the European Database for Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed> and at www.esteticceram.com/downloads.

Basic-UDI-DI: 7630503100007H

For detailed information, please refer the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de

The product properties are based on compliance with and observance of these instructions for use.

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Merz Dental GmbH